

**TECHNICAL REPORT
NATICK/TR-04/011**



AD _____

**IDENTIFICATION OF SELECTED FIBERS
BY PYROLYSIS-GAS CHROMATOGRAPHY-MASS
SPECTROMETRY**

by
Rosa Linda Riel-Bagalawis

April 2004

Final Report
September 2000 – December 2003

Approved for public release; distribution is unlimited

**U.S. Army Research, Development and Engineering Command
Natick Soldier Center
Natick, Massachusetts 01760-5019**

Best Available Copy

20040503 085

DISCLAIMERS

The findings contained in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

Citation of trade names in this report does not constitute an official endorsement or approval of the use of such items.

DESTRUCTION NOTICE

For Classified Documents:

Follow the procedures in DoD 5200.22-M, Industrial Security Manual, Section II-19 or DoD 5200.1-R, Information Security Program Regulation, Chapter IX.

For Unclassified/Limited Distribution Documents:

Destroy by any method that prevents disclosure of contents or reconstruction of the document.

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188													
<p>The public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Department of Defense, Washington Headquarters Services, Directorate for Information Operations and Reports (0704-0188), 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to any penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.</p> <p>PLEASE DO NOT RETURN YOUR FORM TO THE ABOVE ADDRESS.</p>																	
1. REPORT DATE (DD-MM-YYYY) 19 April 2004		2. REPORT TYPE Final		3. DATES COVERED (From - To) September 2000 to December 2003													
4. TITLE AND SUBTITLE IDENTIFICATION OF SELECTED FIBERS BY PYROLYSIS-GAS CHROMATOGRAPHY-MASS SPECTROMETRY				5a. CONTRACT NUMBER													
				5b. GRANT NUMBER													
				5c. PROGRAM ELEMENT NUMBER 423012.12													
6. AUTHOR(S) Rosa Linda Riel-Bagalawis				5d. PROJECT NUMBER OMA													
				5e. TASK NUMBER													
				5f. WORK UNIT NUMBER													
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) U. S. Army Research, Development and Engineering Command Natick Soldier Center ATTN: AMSRD-NSC-IP-E Natick, MA 01760-5019				8. PERFORMING ORGANIZATION REPORT NUMBER NATICK/TR-04/011													
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)													
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)													
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution is unlimited																	
13. SUPPLEMENTARY NOTES																	
14. ABSTRACT A pyrolysis-gas chromatographic method was developed in 2000 to characterize military fibers and fabrics. As the performance and durability of materials relate to their fiber composition, a sensitive and reliable method for the identification of fibers became necessary. The method was further improved by incorporating a mass spectrometer as a detector to ensure accurate component identification from the pyrolyzates. Pyrograms and gas chromatographic analysis are highly reproducible. This analytical method is reliable for the characterization of fibers.																	
15. SUBJECT TERMS <table style="width: 100%; border: none;"> <tr> <td style="width: 33%;">FIBERS</td> <td style="width: 33%;">GAS CHROMATOGRAPHY</td> <td style="width: 33%;">IDENTIFICATION</td> </tr> <tr> <td>FABRICS</td> <td>MASS SPECTROMETRY</td> <td>PYROLYSIS PRODUCTS</td> </tr> <tr> <td>DETECTION</td> <td>ANALYTICAL METHOD</td> <td>POLYMERS</td> </tr> <tr> <td>PYROLYSIS</td> <td>QUALITATIVE ANALYSIS</td> <td>MOLECULAR STRUCTURE</td> </tr> </table>						FIBERS	GAS CHROMATOGRAPHY	IDENTIFICATION	FABRICS	MASS SPECTROMETRY	PYROLYSIS PRODUCTS	DETECTION	ANALYTICAL METHOD	POLYMERS	PYROLYSIS	QUALITATIVE ANALYSIS	MOLECULAR STRUCTURE
FIBERS	GAS CHROMATOGRAPHY	IDENTIFICATION															
FABRICS	MASS SPECTROMETRY	PYROLYSIS PRODUCTS															
DETECTION	ANALYTICAL METHOD	POLYMERS															
PYROLYSIS	QUALITATIVE ANALYSIS	MOLECULAR STRUCTURE															
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT SAR	18. NUMBER OF PAGES 42	19a. NAME OF RESPONSIBLE PERSON Rosa Linda Riel Bagalawis												
a. REPORT U	b. ABSTRACT U	c. THIS PAGE U			19b. TELEPHONE NUMBER (Include area code) (508) 233-4355												

Table of Contents

List of Figures	v
List of Tables	vii
Preface and Acknowledgements	ix
Introduction	1
Materials and Methods	5
Standards	5
Samples	5
Quality Assurance/Quality Control (QA/QC)	5
Pyrolysis-Gas Chromatography (Py-GC)	5
Pyrolysis-Gas Chromatography – Mass Spectrophotometry (Py-GC-MS)	6
Results	7
Conclusions	31
List of References	32

List of Figures

Figure 1.	Pyrogram of cellulose acetate	10
Figure 2.	Mass spectrum of acetic acid	10
Figure 3.	Pyrogram of acrylic (polyacrylonitrile)	11
Figure 4.	Mass spectrum of acrylonitrile	11
Figure 5.	Pyrogram of cotton	12
Figure 6.	Mass spectrum of glycolaldehyde	12
Figure 7.	Pyrogram of Gore-Tex®	13
Figure 8.	Mass spectrum of tetrafluoroethylene	13
Figure 9.	Pyrogram of Kevlar®	14
Figure 10.	Mass spectrum of 1,4-diaminobenzene	14
Figure 11.	Pyrogram of Nomex®	15
Figure 12.	Mass spectrum of 1,3-diaminobenzene	15
Figure 13.	Pyrogram of Nylon 6	16
Figure 14.	Mass spectrum of caprolactam	16
Figure 15.	Pyrogram of nylon 66	17
Figure 16.	Mass spectrum of cyclopentanone	17
Figure 17.	Pyrogram of nylon 12	18
Figure 18.	Mass spectrum of octene	18
Figure 19.	Pyrogram of polyester	19
Figure 20.	Mass spectrum of benzoic acid	19
Figure 21.	Pyrogram of polyethylene	20
Figure 22.	Mass spectrum of decane	20

Figure 23.	Pyrogram of polypropylene	21
Figure 24.	Mass spectrum of 2,4-dimethyl-1-heptene	21
Figure 25.	Pyrogram of rayon	22
Figure 26.	Mass spectrum of glycolaldehyde	22
Figure 27.	Pyrogram of silk	23
Figure 28.	Mass spectrum of phenol	23
Figure 29.	Pyrogram of wool	24
Figure 30.	Mass spectrum of carbon disulfide	24
Figure 31.	Pyrogram of an unknown fabric sample	25
Figure 32.	Mass spectrum of benzoic acid	25
Figure 33.	Pyrogram of an unknown helmet liner sample	26
Figure 34.	Mass spectrum of caprolactam	26
Figure 35.	Pyrogram of an unknown tent material sample after extraction of polyvinyl chloride (PVC)	27
Figure 36.	Mass spectrum of benzoic acid	27
Figure 37.	Pyrogram of an unknown foamy inner layer sample from a headgear	28
Figure 38.	Mass spectrum of 1-decene	28

List of Tables

Table 1.	Molecular Structure of Fibers and Their Major Pyrolysis Products	3
Table 2.	Standard Fibers, Major Pyrolysis Products Identified and their Molecular Weights	7
Table 3.	Samples Identified by Py-GC	29
Table 4.	Samples Identified by Py-GC-MS	30

Preface and Acknowledgements

This analytical work of fibers supports projects in the Engineering, Prototype and Performance Evaluation Team, Army Soldier System Engineering Team and Materials Integrated Team of the Individual Protection Directorate (IPD), U.S. Army Research, Development and Engineering Command (RDECOM), Natick Soldier Center, Natick, Massachusetts. This study also supports projects of the Aerial Delivery Engineering Support Team of the Airdrop/Aerial Delivery Directorate and the 21st Century Fabric Structures Team of the Shelter Directorate belonging to the same Command.

Identification of fibers using a Pyrolysis-Gas Chromatography-Flame Ionization Detector (Py-GC-FID) was conducted by comparing gas chromatographic (GC) profiles of the pyrolysis products of the sample to those of known standards. In 2003, the addition of a Mass Spectrometer (MS) detector to the gas chromatograph (GC) enabled us to confirm the presence of fibers not only from their GC profiles but also by their major pyrolysis products.

The research was conducted at the U. S. Army RDECOM during the period of September 2000 through December 2003 under program element 423012.12 and project OMA.

The author would like to thank Walter Yeomans and Joel Carlson for their help in the configuration of the gas chromatography-mass spectrometry (GC-MS) software after the installation of the pyrolysis apparatus on the new GC-MS.

IDENTIFICATION OF SELECTED FIBERS BY PYROLYSIS- GAS CHROMATOGRAPHY-MASS SPECTROMETRY

Introduction

The U. S. Army utilizes a large number of fibers and fabrics in support of clothing, headgear and tent materials. As the performance and durability of materials relate to their fiber content, a sensitive and reliable method for the identification of fibers became necessary.

Fibers cannot be analyzed by direct gas chromatography and mass spectrometry because they are often blended polymers that are non-volatile. They have to be broken into simpler components that are volatile for analysis.

Pyrolysis is a chemical process whereby thermal energy breaks down large molecules into smaller fragment molecules in the absence of oxygen. Pyrolysis coupled with gas chromatographic (Py-GC) separation of the volatile fragments provides a simple analytical method that can be used to characterize the nature and identity of the fibers and fabrics of military interest. The addition of mass spectrometry to the GC allows one to identify the individual fragments (Py-GC-MS) providing a powerful tool that can be applied to a wide range of chemical problems (Hardin, I., 1996).

The pyrolysis profile is characteristic of a particular sample either in the appearance of unique components or in the relative distribution of the products (Irwin, 1982).

From 2000 - 2002, Natick used Py-GC techniques for the characterization of fibers from clothing, headgear/helmet liners, parachutes, tent materials and other garments. Using Py-GC equipped with a flame ionization detector, we obtained pyrograms (graphs with peaks and retention times) from only a few fibers and based on references, we were able to identify the fibers after matching them to pyrogram profiles in the standards.

For fiber material with less than 5% of additives/finish, sample preparation is not needed for this analysis. Otherwise, extraction of the nonfibrous materials may be necessary as is the case with the analysis of PVC-coated polyester.

In more recent studies, Natick has combined pyrolysis-gas chromatography (Py-GC) with mass spectrometry (MS). For qualitative confirmation purposes, to ensure that the appropriate pyrogram represented a specific fiber, the mass spectrometer was used to identify the major fragments through a mass spectral library search.

Since these processes Py-GC-MS are characteristic of the polymers involved, and are highly reproducible, it becomes relatively straightforward to identify a polymer based on the identity of its pyrolysis products (Wampler and Zawodny, 1999).

Polyolefins such as polyethylene and polypropylene generate many oligomers, creating a characteristic pattern that becomes readily recognizable. Polyethylene is a long, aliphatic polymer; its pyrolytic products are normal hydrocarbons, including alkanes, alkenes and dienes, which produce the pattern of triplet peaks in the pyrogram of polyethylene. The effect of having the methyl branches observed in polypropylene is easily discerned in the pyrogram. The large peak of dimethyl heptene, trimer of propylene is seen in the pyrogram of poly propylene (Wampler and Zawodny, 1999).

By reviewing the structure of a fiber/polymer, the pyrolysis products can be predicted as in the cases of Nomex® and Kevlar®. The only difference between these two polymers is that in Kevlar®, the amino and carboxyl groups are in para positions on the aromatic rings while in Nomex®, the amino and carboxyl groups are in the meta positions of the aromatic rings. In the same way, any fiber with a cellulosic structure can be projected to yield pyrolysis products similar to those of cotton.

Table 1. on the following pages shows molecular structures and the pyrolysis products of selected fibers as cited by references and other products as predicted based on the molecular structures of a few fibers.

Table 1. *Molecular Structure of Fibers and Their Major Pyrolysis Products*

Fiber ID	Molecular Structure	Major Pyrolysis Products
Acetate (Cellulose Acetate)	$\left[\begin{array}{c} \text{OR} \\ \\ -\text{O}-\text{C}_6\text{H}_4-\text{O}-\text{C}(\text{OR})_2-\text{C}_6\text{H}_4-\text{O}- \\ \\ \text{CH}_2\text{OR} \end{array} \right]_n$ <p>(Pouchert, C., 1985)</p> <p>R=CCH₃</p>	Acetic anhydride; acetic acid; glycolaldehyde (Hardin, I. & Wang, X. 1989)
Acrylic (Polyacrylonitrile)	$\left[-\text{CHCH}_2- \right]_n$ <p>(Pouchert, C., 1985)</p> <p>C≡N</p>	Acetonitrile; acrylonitrile or propenenitrile (Wampler, T. & Zawodny, C., 1999)
Cotton	$\left[\begin{array}{c} \text{OH} \\ \\ -\text{O}-\text{C}_6\text{H}_4-\text{O}-\text{C}(\text{OH})_2-\text{C}_6\text{H}_4-\text{O}- \\ \\ \text{CH}_2\text{OH} \end{array} \right]_n$ <p>(Segal, L. & Wakelyn, P. J., 1985)</p>	Carbon dioxide; oxirane; furan; acetone; glycolaldehyde; acetic acid; 1-hydroxy 2-propanone; acetic anhydride; furfuraldehyde (Hardin, I., 1996)
Gore-Tex® (PTFE)	$\left[\begin{array}{c} \text{F} \quad \text{F} \\ \quad \\ -\text{C}-\text{C}- \\ \quad \\ \text{F} \quad \text{F} \end{array} \right]_n$ <p>(Stecher, P., Finkel, M., Siegmund, O. & Szafranski, B., 1960)</p>	Tetrafluoroethylene (Wampler, T. & Zawodny, C., 1999)
Kevlar®	$\left[\begin{array}{c} \text{NH} \quad \text{CO} \\ \quad \\ -\text{C}_6\text{H}_4-\text{NH}-\text{CO}-\text{C}_6\text{H}_4-\text{NH}-\text{CO}- \\ \quad \\ \text{NH} \quad \text{CO} \end{array} \right]_n$ <p>(Yang, H., 1993)</p>	Benzene; aniline; benzonitrile; amino compound, carbonitriles (Wampler, T., 1995); 1,4-benzenediamine (projected)
Nomex®	$\left[\begin{array}{c} \text{O} \quad \text{O} \\ \quad \\ -\text{C}_6\text{H}_4-\text{N}-\text{C}-\text{C}_6\text{H}_4-\text{N}-\text{C}- \\ \quad \\ \text{H} \quad \text{H} \end{array} \right]_n$ <p>(Yang, H., 1993)</p>	Benzene; toluene; aniline; benzonitrile; amino compound; carbonitriles (Wampler, T., 1995); 1,3-benzenediamine (projected)

Table 1. *Molecular Structure of Fibers and Their Major Pyrolysis Products (Cont'd)*

Fiber ID	Molecular Structure	Major Pyrolysis Products
Nylon 6	$\begin{array}{c} \text{O} \\ \parallel \\ [-\text{NHCH}_2(\text{CH}_2)_4\text{CH}_2\text{C}-]_n \\ \text{(Pouchert, C., 1985)} \\ \text{-NH} \\ \\ \text{CH}_2 \\ \\ (\text{CH}_2)_4 \\ \\ \text{CH}_2 \\ \\ \text{O} \\ \parallel \\ [-\text{NH}-\text{CCH}_2(\text{CH}_2)_4\text{CH}_2\text{C}-]_n \\ \text{(Pouchert, C., 1985)} \end{array}$	Caprolactam (CDS Applications Lab Personnel, 1995)
Nylon 66	$\begin{array}{c} \text{O} \\ \parallel \\ [-\text{NHCH}_2(\text{CH}_2)_4\text{CH}_2\text{C}-]_n \\ \text{(Pouchert, C., 1985)} \\ \text{O} \\ \parallel \\ [-\text{NH}-\text{CCH}_2(\text{CH}_2)_3\text{CH}_2\text{C}-]_n \\ \text{(Pouchert, C., 1985)} \end{array}$	Cyclopentanone; hexamethylene diamine (Wampler, T., 1995)
Nylon 12	$\begin{array}{c} \text{O} \\ \parallel \\ [-\text{NHCH}_2(\text{CH}_2)_9\text{CH}_2\text{C}-]_n \\ \text{(Pouchert, C., 1985)} \end{array}$	Amides; hydrocarbons; mononitriles (CDS Application Lab Personnel, 1995)
Polyester	$\begin{array}{c} \text{O} \\ \parallel \\ [-\text{O}(\text{CH}_2)_4\text{OC}-\text{C}_6\text{H}_4-\text{C}-]_n \\ \text{(McIntyre, J. E., 1985)} \end{array}$	Benzene; benzoic acid; biphenyl; vinyl terephthalate (Hardin, I., 1996); benzene, vinyl benzoate, benzoic acid, biphenyl (Wampler, T., 1995)
Polyethylene	$[-\text{CH}_2\text{CH}_2-]_n$ <p>(Pouchert, C., 1985)</p>	1,9-decadiene; 1-decene; 1-decane (C10 up to C39) (Wampler, T., 1995)
Polypropylene	$[-\text{CH}_2\text{CH}(\text{CH}_3)-]_n$ <p>(Pouchert, C., 1985)</p>	Methyl branched hydrocarbons: trimer, tetramer, etc. Dimethylheptene (Wampler, T. & Zawodney, C., 1999)
Rayon	$\left[\begin{array}{c} \text{CH}_2\text{OH} \quad \text{CH}_2\text{OH} \\ \quad \\ \text{C}_6\text{H}_4 \quad \text{C}_6\text{H}_4 \\ \quad \\ \text{OH} \quad \text{OH} \end{array} \right]_n$ <p>(Dyer, J. & Daul, G., 1985)</p>	Carbon dioxide; acetic acid; glycolaldehyde; 1-hydroxy-2-propanone; furaldehyde; 1,2-cyclopentanedione (projected)
Silk	$\begin{array}{c} \text{O} \\ \parallel \\ \text{H}_2\text{N}-\text{CH}_2-\text{C}-\text{OH} \\ \\ \text{Glycine} \\ \text{Plus other amino acids (Pouchert, C., 1985)} \end{array}$	Benzene, toluene, picoline (Irwin, 1982) Carbon dioxide; aromatic compound; phenol (CDS Application Lab Personnel, 1995)
Wool	$\begin{array}{c} \text{O} \\ \parallel \\ \text{H}_2\text{N}-\text{CH}_2-\text{C}-\text{OH} \\ \\ \text{Glycine} \\ \text{Plus other amino acids (Pouchert, C., 1985)} \end{array}$	Carbon dioxide; methane; hydrogen sulfide; Carbon disulfide (Hardin, I., 1996)

MATERIALS AND METHODS

Standards

Standard fibers of acetate, acrylic, cotton, Gore-Tex®, Kevlar®, Nomex®, nylon 6, nylon 66, nylon 12, polyester, polyethylene, polypropylene, rayon, silk and wool were obtained from our collection here at Natick. Approximately a microgram of each standard was required for each study.

Samples

A microgram of each sample was required and the number of samples varied from 2 to 23.

Quality Assurance/Quality Control (QA/QC)

Prior to any analysis, the pyrolysis chamber was auto-cleaned by setting it to 1200°C for 20 seconds as recommended by the equipment manufacturer. An instrument blank was run on the GC-MS to ensure that the instrument was clear of all interferences.

Laboratory blanks, in which the pyrolyzer was operated with no sample present, was used to ensure that each analysis was free of contamination. Any blank sample showing a positive concentration resulted in the qualification of the analytical sample. In these studies, no contaminated blanks were observed.

Standard duplicates were also used to ensure the reproducibility of the analytical method. In the event that the pyrogram failed to yield reproducible results, the analytical instrument was diagnosed and repaired or sent out for preventive maintenance or repair.

Auto-cleaning of the pyrolysis chamber was conducted after every sample. Standards were run every 12 hours of operation.

Pyrolysis-Gas Chromatography (Py-GC)

The first Py-GC technique was developed using a Model AS 2500 Pyrolyzer equipped with an auto sampler manufactured by CDS Analytical connected to a Hewlett-Packard Model 5890A Gas Chromatograph. The pyrolyzer chamber/oven and the auto sampler and gas chromatograph were all computer controlled. About one microgram of a sample was placed into

20 seconds. The pyrolyzer interface temperature that transferred the pyrolyzates to the GC was set at 300°C. The probe was automatically cleaned after every run at 1200°C for 20 seconds.

The carrier gas Helium flowing 5 mLs per minute swept the volatile through the splitless inlet of the gas chromatograph at a temperature of 300°C and into a Hewlett-Packard HP-5 column of fused silica, 30-meters long with 0.53 mm I.D. coated with 5%-Phenyl-methylpolysiloxane up to the detector. The oven temperature was initially maintained at 40°C for one minute and then it was temperature programmed at the rate of 6° per minute up to a final temperature of 300°C for 15 minutes. Hydrogen and compressed air were used for the flame ionization detector that was maintained at 275°C.

Peak heights, areas and retention times for the standards and samples were measured and sequentially saved by a computer. The pyrograms of the samples were compared with those of the standards. The characterization of the sample fibers was conducted by comparing the sample pyrograms with those of the standards.

Pyrolysis-Gas Chromatography-Mass Spectrometry(Py-GC-MS)

In recent studies, the Model AS 2500 Pyrolyzer manufactured by CDS Analytical was combined with a Gas Chromatograph Model 6890N- Mass Spectrometer Model 5873N in order to enhance the selectivity and sensitivity of the analytical method. The column was changed from HP-5 to HP-5ms column of bonded 5%-Phenyl-methylpolysiloxane, 30-meters long and 0.25 mm I.D. The carrier gas was helium but hydrogen and air were no longer required with the elimination of the Flame Ionization Detector from the method. The chromatographic conditions described above were used on the GC-MS except the flow rate of helium was reduced to 1 mL per minute to accommodate the requirements of the mass spectrometer. The detector temperature was 280°C. The inlet flow split ratio was 50:1.

Results

The pyrolysis products of the standard fibers were identified and their molecular weights were obtained searching the NIST Mass Spectral Library (Stein, Levitsky, Fateev, Tchekhovski & Mallard, 2000). See Table 2. below.

Table 2. *Standard Fibers, Major Pyrolysis Products Identified and their Molecular Weights*

Fiber ID	Major Pyrolysis Products Identified	Mol. Wt.
Acetate (Cellulose acetate)	Acetic acid	60
	Methyl acetate	74
Acrylic (Polyacrylonitrile)	Propenenitrile (acrylonitrile)	53
Cotton (cellulose-based)	Carbon dioxide	44
	Glycolaldehyde	60
	1-hydroxy-2-propanone	74
	Furfuraldehyde	96
	2-hydroxy-2-cyclopentene-1-one	112
Gore-tex® (Teflon)	Tetrafluoroethene or tetrafluoroethylene	100
Kevlar® (p-phenylene terephthalamide)	Aniline	93
	Benzonitrile	103
	1,4-benzenediamine	108
Nomex® (m-phenylene isophthalamide)	Aniline	93
	Benzonitrile	103
	1,3-benzenediamine	108
Nylon 6	Caprolactam	113
Nylon 66	Cyclopentanone	84

Table 2. *Standard Fibers, Major Pyrolysis Products Identified and their Molecular Weights (Cont')*

Fiber ID	Major Pyrolysis Products Identified	Mol. Wt.
Nylon 12	2-propenenitrile	53
	Nonene	126
	Decanenitrile	153
Polyester (PET)	Benzene	78
	Benzoic acid	122
	Biphenyl	154
Polyethylene	1,9-decadiene	138
	Decene	140
	Decane	142
	Triplets of -diene, -olefin & n-alkane, C10-C39	
Polypropylene	2-methyl-1-pentene	84
	2,4-dimethyl-1-heptene	126
Rayon (cellulose-based)	Carbon dioxide	44
	Acetic acid	60
	Glycolaldehyde	60
	1-hydroxy-2-propanone	74
	Furaldehyde	96
Silk (protein)	Carbon dioxide	44
	Acetyl hydrazide	74
	Ammonium carbamate	78
	Phenol	94
	Methyl phenol	108
Wool (sulfur protein)	Carbon disulfide	76
	Ammonium carbamate	78
	Toluene	92
	Phenol	94
	4-methyl phenol	108

On the pages that follow, Figures 1 to 30 show pyrograms of selected standard fibers and mass spectra of their major pyrolysis products. By comparing the pyrogram and the major pyrolysis products of a sample with the pyrogram and the major pyrolysis products of the known standard, the unknown fabric can be identified. The following conclusions can be made using the data collected. Figures 31 and 32 on the following pages show a pyrogram of an unknown fabric material sample and a mass spectrum of benzoic acid, respectively. The pyrogram of the sample and the major pyrolysis product, benzoic acid are comparable to that of polyester. Figures 33 and 34 show a pyrogram of an unknown helmet liner sample and mass spectrum of caprolactam, the major pyrolysis product, respectively. The pyrogram and the major pyrolysis products are comparable to those of nylon 6. Figure 35 shows a pyrogram of a tent material sample after extraction of polyvinyl chloride. It is similar to that of polyester. Figure 36 is a mass spectrum of benzoic acid, a major pyrolysis product of polyester; therefore the tent material sample is a mixture of polyester and polyvinyl chloride. Figure 37 shows a pyrogram of an unknown foamy inner layer of a headgear sample. The profile is similar to that of polyethylene. One of the major pyrolysis products is 1-decene as shown by Figure 38. The other pyrolysis products were 1,9-decadiene and decane, therefore it is polyethylene.




Figure 1. Pyrogram of cellulose acetate




Figure 2. Mass spectrum of acetic acid




Figure 3. Pyrogram of polyacrylonitrile




Figure 4. Mass spectrum of acrylonitrile




Figure 5. Pyrogram of cotton




Figure 6. Mass spectrum of glycolaldehyde




Figure 7. Pyrogram of Gore-tex




Figure 8. Mass spectrum of tetrafluoroethylene




Figure 9. Pyrogram of Kevlar®




Figure 10. Mass spectrum of 1,4-diaminobenzene




Figure 11. Pyrogram of Nomex®




Figure 12. Mass spectrum of 1,3-diaminobenzene




Figure 13. Pyrogram of nylon 6




Figure 14. Mass spectrum of caprolactam




Figure 15. Pyrogram of nylon 66




Figure 16. Mass spectrum of cyclopentanone




Figure 17. Pyrogram of nylon 12




Figure 18. Mass spectrum of 2-propenenitrile




Figure 19. Pyrogram of polyester




Figure 20. Mass spectrum of benzoic acid




Figure 21. Pyrogram of polyethylene




Figure 22. Mass spectrum of decene




Figure 23. Pyrogram of polypropylene




Figure 24. Mass spectrum of 2,4-dimethyl-1-heptene




Figure 25. Pyrogram of rayon




Figure 26. Mass spectrum of glycolaldehyde




Figure 27. Pyrogram of silk




Figure 28. Mass spectrum of phenol




Figure 29. Pyrogram of wool




Figure 30. Mass spectrum of carbon disulfide

Figures 31 to 38 show pyrograms of samples and mass spectra of their major pyrolysis products.




Figure 31. Pyrogram of an unknown fabric sample

The pyrogram above in Figure 31 and the mass spectrum in Figure 32 confirms the presence of polyester.




Figure 32. Mass spectrum of benzoic acid




Figure 33. Pyrogram of an unknown helmet liner sample

The pyrogram above in Figure 33 and the mass spectrum in Figure 34 confirm the presence of nylon 6.




Figure 34. Mass Spectrum of Caprolactam




Figure 35. Pyrogram of an unknown tent material sample after extraction of the polyvinyl chloride (PVC)

The pyrogram above in Figure 35 and the mass spectrum in Figure 36 confirm the presence of polyester.




Figure 36. Mass spectrum of benzoic acid




Figure 37. Pyrogram of an unknown foamy inner layer sample from a headgear

The pyrogram above in Figure 37 and the mass spectrum in Figure 38 confirm the presence of polyethylene.




Figure 38. Mass spectrum of 1-decene

Not all mass spectra of major pyrolysis products are shown here because there are just too many. However, they are summarized on Tables 2 as matched with the mass spectral library showing corresponding molecular weights for the standards and on Table 3 below and Table 4 on next page for the samples.

Table 3. *Samples Identified by Py-GC*

Sample Number	Matched Pyrogram Profile	Major Pyrolysis Product*	Fiber Identified
8415-01-444-1169	Nylon 66	Cyclopentanone	Nylon 66
8415-01-444-1200	Nylon 66	Cyclopentanone	Nylon 66
8415-01-444-2325	Nylon 66	Cyclopentanone	Nylon 66
8415-01-444-2308	Nylon 66	Cyclopentanone	Nylon 66
BL15019832 4HBF	Nylon 66	Cyclopentanone	Nylon 66
BL15019832 4H	Nylon 66	Cyclopentanone	Nylon 66
090597103TOLAH	Nylon 66	Cyclopentanone	Mostly Nylon 66
090597103TOWHBF	Nylon 66	Cyclopentanone	Mostly Nylon 66
090597103TOWH	Nylon 66	Cyclopentanone	Mostly Nylon 66
BL261197232TOPBH	Nylon 66	Cyclopentanone	Mostly Nylon 66
BL261197232TOBHW	Nylon 66	Cyclopentanone	Mostly Nylon 66
BL261197232TOWEB	Nylon 6	Caprolactam	Nylon 6
BL261197232TOW	Nylon 66	Cyclopentanone	Mostly Nylon 66
B212059729OCP1J	Nylon 66	Cyclopentanone	Mostly Nylon 66
1JHBN2A2AUTOCLAV	Nylon 66	Cyclopentanone, etc.	Nylon 66 +
P1P11P1030991HBP	Nylon 66	Cyclopentanone, etc.	Nylon 66 +
P1P11P1030991HBR	Nylon 66	Cyclopentanone, etc.	Nylon 66 +
P1P11P1030991HRP	Nylon 66	Cyclopentanone, etc.	Nylon 66 +
P1P11P1030991SHP	Nylon 66	Cyclopentanone, etc.	Nylon 66 +
016F1601	Rayon	Not identified	Rayon
011F0111	Polyester	Benzoic Acid	Polyester
015F0101	Nylon 66	Cyclopentanone	Nylon 66
MEMB0019	Polyethylene	3 Hydrocarbons	Polyethylene
PM000001	Polyester	Benzoic Acid	Polyester

* The major pyrolysis product was identified by position and height on the pyrogram as identified on references. The fiber was identified by pyrogram profile comparison of sample to those of standards.

Table 4. *Samples Identified by Py-GC-MS*

Sample Number	Matched Pyrogram Profile	Major Pyrolysis Product(s) & Mol. Wt.	Fiber Identified
04220302	Acrylic	Acrylonitrile - 53	Acrylic
04220303	Nylon 66	Cyclopentanone - 84	Nylon 66
04220305	Polyester	Benzene - 78 Benzoic acid - 122	Polyester
070803	Nylon 6	Caprolactam - 113	Nylon 6
0709031	Polyester	Benzoic acid - 122 Biphenyl - 154	Polyester
0709032	Polyester	Benzoic Acid - 122 Biphenyl - 154	Polyester
0709033	Nylon 6	Caprolactam - 113	Nylon 6
07160301	Polyethylene	1-Decene - 140 Decane - 142	Polyethylene
07160302	Polyethylene	1,9-Decadiene - 138 1-Decene - 140 Decane - 142	Polyethylene
07160303	Nylon 66	Cyclopentanone - 84	Nylon 66
07160304	Nylon 6	Caprolactam - 113	Nylon 6
9250301	Polyester	Benzene - 78 Benzoic Acid - 122 Biphenyl - 154	Polyester
9250302	Polyester	Benzene - 78 Benzoic acid - 122	Polyester

Conclusions

The addition of the mass spectrometer to the GC has expanded our capability to identify unknown fibers. When MS is used in conjunction with Py-GC, more conclusive results can be obtained in identifying unknown materials.

A pyrogram is considered a "fingerprint" of a certain fiber. It is very much reproducible as long as the conditions are kept constant. Additionally, when a major peak(s) from the pyrolysis of a sample is/are identified by mass spectrometry to be the same as that/those of the standard, there is no more doubt that they are the same fibers.

Py-GC-MS is an excellent tool for the characterization of fibers used in military materials and has allowed a more thorough analysis of fibers/fabrics used in military applications.

This document reports research undertaken at the U.S. Army Research, Development and Engineering Command, Natick Soldier Center, Natick, MA, and has been assigned No. NATICK/TR-04/011 in a series of reports approved for publication.

References

- Bergen, V. and Mauersberger, H. R. American Wool Handbook, Second Edition, 174. Textile Book Publishers, Inc. New York, New York. 1948.
- CDS Application Lab Personnel. "A Comparison of Nylons by Pyrolysis-GC," 1. CDS Analytical, Inc. Oxford, Pennsylvania. 1995.
- CDS Application Lab Personnel. "Pyrolysis-GC of Natural Fiber," 1. CDS Analytical, Inc. Oxford, Pennsylvania. 1995.
- Dyer, J. and Daul, G. C. "Rayon Fibers." Lewin, M. and Pearce, E. M. Handbook of Fiber Science and Technology: Fiber Chemistry, Volume IV, 909. Marcel Dekker, Inc. New York, New York. 1985.
- Hardin, I. R. "Investigation of Textiles by Analytical Pyrolysis." Raheel, M. Modern Textile Characterization Methods, 175. Marcel Dekker, Inc. New York, New York. 1996.
- Hardin, I. & Wang, X. "The Use of Pyrolysis-Gas Chromatography in Textiles as an Identification Method." Textile Colorist and Chemist, Volume 21, Number 1, 29. American Association of Textile Chemists and Colorists. Research Triangle Park, North Carolina. 1989.
- Irwin, W. J. Analytical Pyrolysis, A Comprehensive Guide, Volume 22, 334. Marcel Dekker, Inc. New York, New York. 1982.
- McINTYRE, J. E. "Polyester Fibers." Lewin, M. and Pearce, E. M. Handbook of Fiber Science and Technology: Fiber Chemistry, Volume IV, 1. Marcel Dekker, Inc. New York, New York. 1985.
- Pouchert, C. Aldrich Library of FT-IR Spectra, Edition I, Volume 2, 1157. Aldrich Chemical Co. Inc. Milwaukee, Wisconsin. 1985.
- Robson, R. M. "Silk." Lewin, M and Pearce, E. M. Handbook of Fiber Science and Technology: Fiber Chemistry, Volume IV, 647. Marcel Dekker, Inc. New York, New York. 1985.
- Segal, L. & Wakelyn, P. J. "Cotton Fibers." Lewin, M & Pearce E. M. Handbook of Fiber Science and Technology: Fiber Chemistry, Volume IV, 809. Marcel Dekker, Inc. New York, New York. 1985.
- Stecher, P., Finkel, M., Siegmund, O. & Szafranski, B. The Merk Index of Chemicals and Drugs, 7th Edition, 834. Merk Index & Company. Rahway, New Jersey. 1960.

- Stein, Levitsky, Fateev, Tchekhovski, and Mallard. NIST Mass Spectral Libraries. 101. Agilent Technologies Technical Publication Department. Palo Alto, California. 2000
- Wampler, T. Applied Pyrolysis Handbook, 13. Marcel Dekker, Inc. New York, New York. 1995.
- Wampler, T. & Zawodny, C. "*Analysis of polymer packaging products using pyrolysis-gas chromatography-mass spectrometry.*" American Laboratory News Edition, Volume 31, Number 19, Application Note, 30. American Laboratory. Sheldon, Connecticut. 1999.
- Yang, H. H. "*Nomex® Aramid Fiber.*" Lewin, M. & Preston, J. Handbook of Fiber Science and Technology: High Technology Fibers, Volume III, Part C, 98. Marcel Dekker, Inc. New York, New York. 1993.